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mark	
eegDatv	
headerLength.....	length of header
fileLengthBytes	length of rawdata in bytes
fileLengthSec	length of rawdata in seconds

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```
keyNr ..... S/W key #
gf ..... generalAmpFeatures
spare
```

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```

amplifierID
subAmpID
maxNchAmp ..... max number of channels
digiChannel..... digital/photo channel
epVersion..... ep-data version of data
differentAdcTypes ..... # of different ADCs within amplifier channels
differentGain ..... # of different Gains within amplifier channels
differentHP ..... # of different HP within amplifier channels
differentTP ..... # of different TP within amplifier channels
spare

```

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```

nrChan..... number of channels in data file
source..... 1...ONLINE created
             2...OFFLINE created
nChAd..... # AD-channels = NchAmp + NchDig_activ
nChAmp..... # Amplifier channels
nChDig_activ ..... # Digital channels
nChDsk; ..... # stored channels: NchAmp + NchDig + 1 (Marker)
nHz ..... samplefrequency
nPps ..... points a second
nrInvasiv ..... invasive electrodes included
dummy ..... spare
info..... amplifier info

```

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hp..... HP [Hz]
 tp TP [Hz]
 chanGain Gain
 notch..... 0 or 1 (for 50/60 Hz)
 res1..... reserve

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act[6]..... name of active electrode
 ref[6] name of reference electrode
 res0..... spare
 mikroVperBin ADC: $\mu\text{V}/\text{Bin}$
 chanTyp..... mono/bipol/etc
 nrBits ADC: 12 or 16 Bit
 pos3Dxyz position of electrodes (from digitizer, etc.)
 res1..... spare

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